

PRAGATHI MEDICAL SYSTEMS

Premier Resuscitation Unit - PRU2101

FEATURES AND SPECIFICATIONS

- Microprocessor Based Control
- Dual Probe: Skin/Air
- 4 Modes of Operation: Baby/Air/Manual/Pre-Warm
- Temperature Range: 20 c to 40 c with 0.1 c Resolution
- Manual Mode With Adjustable Heater Output
- Heater Output: 0% to 100% in increments of 10%
- 3 Dedicated Displays For Set Temperature, Baby/Air Temperature, and Heater Output
- Ceramic Infrared Heater Enclosed In Stainless Steel Parabolic Reflector
- Three Level Safety Cutoff, "Baby Fully Safe" Principle
- Sensor - Heater - Power Failure Alarm
- High/Low Temperature From Set Temperature
- Baby Feed Timer (configurable)
- Feather Touch Keys
- Baby Cradle Head Up/Down with Mattress
- 3 Sides Collapsible Acrylic Walls
- LED Observation Bulb
- Power Requirement: 220V/50Hz
- Rotatable Heater Box with Topside Phototherapy
- Monitor Tray and IV Stand
- Rigid Castor Wheels - 2+2 (2 Brake)
- Body in Mild Steel/Powder Coated
- Down Side Phototherapy
- X-Ray Cassette Holder
- Dual Lower Drawers
- Dual Oxygen Cylinder Holder
- Patient Monitor - Contec CMS8000
- Pulse Oximeter - Bionet Oxy9wave
- Syringe Pump - Hawkmed HK- 400
- Infusion Pump - Hawkmed HK - 100 II
- Inbuilt Slow Suction
- Inbuilt Oxygen Flow Meter
- Oxygen Hood

FRONT PANEL



E-mail: sales@pragathineo-natal.com

Phone: +91-9844111792

+91-9449764764





PREMIUM WARMER

PWU2102

Micro Processor based controller,
With dual probes: Air , Skin.
Three Modes of operation: Baby/Air/Manual
Three Level Safety Cutoff
"Baby Fully Safe Principle"
Fibre Based Body

Features & Specs

- Inbuilt Weighing Scale (Optional) ✓
- X-ray cassette Holder (Optional) ✓
- Oval Shaped Baby Cradle ✓

More Information

CONTACT US

+91-9844111792

sales@pragathineo-natal.com



www.pragathineo-natal.com



Premium Warmer PWU2102

- 1 Microprocessor based controller
- 2 Dual Probes: Air / Skin
- 3 Three modes of operation: Air / Baby / Manual
- 4 Three level safety cutoff "Baby fully safe principle"
- 5 Ceramic IR heater enclosed in SS parabolic reflector
- 6 Fibre Body
Oval shaped baby cradle
- 7 X-ray cassette holder (optional)
Oxygen cylinder holder (optional)





INFANT RADIANT WARMER

SBW2103

FEATURES AND SPECS

- Microprocessor based control
- Dual Probe: Air/Skin
- 3 modes of operation: Air/Baby/Manual
- Temperature Range: 20°C to 40°C
- Resolution: 0.1°C
- Manual Mode with adjustable heater output
- Heater output: 20% to 100% in increments of 20%
- Dedicated Displays for Set and Baby/Air temperature
- Ceramic IR heater enclosed in SS Parabolic reflector
- Feather touch keys
- LED observation bulb
- Power Requirement: 230v AC/50 Hz
- Alarms
 - Power Failure
 - Sensor Failure
 - Heater Failure
 - High Temperature
 - Low Temperature

PHYSICAL CHARACTERISTICS

- Baby Cradle Head Up - Down
- 3 side collapsible acrylic side walls
- Rotatable Heater Box
- Monitor tray and IV Stand
- Rigid Castor Wheels 2+2 (Brakes)
- Body in Mild Steel - Powder Coated



+91-9844111792



sales@paragathineo-natal.com



Since 1996



Features:

- Microprocessor based controller
- Dual Probe: Air / Skin
- 3 Modes of operation: Air/Baby/Manual
- Top-surface hi-intensity LED phototherapy 12+12 LED
- Under-surface hi-intensity LED phototherapy 15 LED
- LED observation light

Infant Warmer with Double Surface LED Phototherapy SBW2105

Ph: +91-9844111792
sales@pragathineo-natal.com

www.pragathineo-natal.com



INFANT RADIANT WARMER 2106, 2106A

Micro Processor based controller,
With dual probes: Air , Skin.

Three Modes of operation:

Baby/Air/Manual

Three Level Safety Cutoff
"Baby Fully Safe Principle"

Features & Specs

Monitor Tray and IV Stand ✓

Detachable Baby Trolley ✓

Rotatable Heater Box ✓

More Information

CONTACT US

+91-9844111792

sales@pragathineo-natal.com



PRAGATHI MEDICAL SYSTEMS

Stand Warmer SBW 2106

SINCE 1996
25
YEARS
OF
DEDICATION
www.pragathineo-natal.com



FEATURES AND SPECIFICATIONS

- Microprocessor Based Control
- Dual Probe: Skin/Air
- 3 Modes of Operation: Baby/Air/Manual
- Temperature Range: 20 c to 40 c with 0.1 c Resolution
- Manual Mode With Adjustable Heater Output
- Heater Output: 20% to 100% in increments of 20%
- 2 Dedicated Displays For Set Temperature and Baby/Air Temperature
- Ceramic Infrared Heater Enclosed In Stainless Steel Parabolic Reflector
- Three Level Safety Cutoff, "Baby Fully Safe" Principle
- Sensor - Heater - Power Failure Alarm
- High/Low Temperature From Set Temperature
- Feather Touch Keys
- LED Observation Bulb
- Power Requirement: 220V/50Hz
- Rotatable Heater Box
- Monitor Tray and IV Stand
- Rigid Castor Wheels - 2+2 (2 Brake)
- Body in Mild Steel/Powder Coated

E-mail: sales@pragathineo-natal.com

Phone: +91-9844111792

+91-9449764764





Simple Baby Warmer - SBW 2108

FEATURES AND SPECIFICATIONS

- Microprocessor Based Control
- Dual Probe: Skin/Air
- 3 Modes of Operation: Baby/Air/Manual
- Temperature Range: 20 c to 40 c with 0.1 c Resolution
- Manual Mode With Adjustable Heater Output
- Heater Output: 20% to 100% in increments of 20%
- 2 Dedicated Displays For Set Temperature and Baby/Air Temperature
- Ceramic Infrared Heater Enclosed In Stainless Steel Parabolic Reflector
- Three Level Safety Cutoff, "Baby Fully Safe" Principle
- Sensor - Heater - Power Failure Alarm
- High/Low Temperature From Set Temperature
- Feather Touch Keys
- Baby Cradle Head Up/Down with Mattress
- 3 Sides Collapsible Acrylic Walls
- LED Observation Bulb
- Power Requirement: 220V/50Hz
- Rotatable Heater Box
- Monitor Tray and IV Stand
- Rigid Castor Wheels - 2+2 (2 Brake)
- Body in Mild Steel/Powder Coated



13485:2016



Since 1996

Deluxe Baby Warmer

DBW2109

Features & Specs

- ✓ Microprocessor based controller
- ✓ Dual probe: Air / Skin
- ✓ Three modes of operation: Baby / Air/ Manual
- ✓ 3 Dedicated displays for Air/Baby temperature
Set temperature
Heater output
- ✓ Ceramic IR Heater enclosed in SS parabolic reflector
- ✓ Three level safety cutoff "Baby Fully Safe Principle"
- ✓ Three side collapsible acrylic side walls (Hinge Mechanism - Rotatable)
- ✓ Rotatable heater box
- ✓ Castors with brakes (2+2)
- ✓ Downside storage tray
- ✓ Monitor tray and IVC holder
- ✓ Power req. - 220 V 50 Hz



www.pragathineo-natal.com

+91-9844111792
 sales@pragathineo-natal.com

PRAGATHI MEDICAL SYSTEMS

Deluxe Baby Warmer with Downside LED Phototherapy DBW 2110

FEATURES AND SPECIFICATIONS

- Microprocessor Based Control
- Dual Probe: Skin/Air
- 3 Modes of Operation: Baby/Air/Manual
- Temperature Range: 20 c to 40 c with 0.1 c Resolution
- Manual Mode With Adjustable Heater Output
- Heater Output: 20% to 100% in increments of 20%
- Dedicated Display For Set Temperature, Baby/Air Temperature and Heater output
- Ceramic Infrared Heater Enclosed In Stainless Steel Parabolic Reflector
- Three Level Safety Cutoff, "Baby Fully Safe" Principle
- Sensor - Heater - Power Failure Alarm
- High/Low Temperature From Set Temperature
- Feather Touch Keys
- Baby Cradle Head Up/Down with Mattress
- 3 Sides Collapsible Acrylic Walls
- LED Observation Bulb
- Rotatable Heater Box
- Monitor Tray and IV Stand
- Rigid Castor Wheels - 2+2 (2 Brake)
- Power Requirement: 220V 50Hz
- Body in Mild Steel/Powder Coated
- Down Side Phototherapy with 15 Bright Blue LEDs
- X-Ray Cassette Holder (Optional)
- Lower Drawers (Optional)
- Oxygen Cylinder Holder (Optional)



E-mail: sales@pragathineo-natal.com
Phone: +91-9844111792
+91-9449764764



13485:2016

PRAGATHI MEDICAL SYSTEMS



Deluxe Baby Warmer with Doubleside LED Phototherapy DBW 2111

FEATURES AND SPECIFICATIONS

- Microprocessor Based Control
- Dual Probe: Skin/Air
- 3 Modes of Operation: Baby/Air/Manual
- Temperature Range: 20c to 40c with 0.1c Resolution
- Manual Mode With Adjustable Heater Output
- Heater Output: 20% to 100% in increments of 20%
- Dedicated Display For Set Temperature, Baby/Air Temperature and Heater output
- Ceramic Infrared Heater Enclosed In Stainless Steel Parabolic Reflector
- Three Level Safety Cutoff, "Baby Fully Safe" Principle
- Sensor - Heater - Power Failure Alarm
- High/Low Temperature From Set Temperature
- Feather Touch Keys
- Baby Cradle Head Up/Down with Mattress
- 3 Sides Collapsible Acrylic Walls
- LED Observation Bulb
- Rotatable Heater Box
- Monitor Tray and IV Stand
- Rigid Castor Wheels - 2+2 (2 Brake)
- Power Requirement: 220V 50Hz
- Body in Mild Steel/Powder Coated
- Down Side Phototherapy with 15 Bright Blue LEDs
- Top Side Phototherapy with 12 +12 Bright Blue LEDs
- X-Ray Cassette Holder (Optional)
- Lower Drawers (Optional)
- Oxygen Cylinder Holder (Optional)

E-mail: sales@pragathineo-natal.com
Phone: +91-9844111792
+91-9449764764



13485:2016

PRAGATHI MEDICAL SYSTEMS

Deluxe Baby Warmer III DBW 2109B

FEATURES AND SPECIFICATIONS

- Microprocessor Based Control
- Dual Probe: Skin/Air
- 3 Modes of Operation: Baby/Air/Manual
- Temperature Range: 20 c to 40 c with 0.1 c Resolution
- Manual Mode With Adjustable Heater Output
- Heater Output: 20% to 100% in increments of 20%
- Dedicated Display For Set Temperature, Baby/Air Temperature and Heater output
- Quarts / Ceramic IR Heater Enclosed In Stainless Steel Parabolic Reflector
- Three Level Safety Cutoff, "Baby Fully Safe" Principle
- Sensor - Heater - Power Failure Alarm
- High/Low Temperature From Set Temperature
- Feather Touch Keys
- Baby Cradle Head Up/Down with Mattress
- 3 Sides Collapsible Acrylic Walls
- LED Observation Bulb
- Rotatable Heater Box
- Monitor Tray and IV Stand
- Rigid Castor Wheels - 2+2 (2 Brake)
- Power Requirement: 220V 50Hz
- Body in Mild Steel/Powder Coated
- Top surface LED Phototherapy
- X-Ray Cassette Holder
- Lower Drawers
- Oxygen Cylinder Holder



E-mail: sales@pragathineo-natal.com
Phone: +91-9844111792
+91-9449764764

